


INTERACTIVE EDITION

Drivers & Tools Catalog 2024/2025




HAND DRIVERS

FOR PROSTHETICS

	HEX (mm)	L (mm)	CODE
	1.25mm	7	HHD1207
		10	HHD1210
		15	HHD1215
		30	HHD1230

HEX DRIVERS


FOR PROSTHETICS

	HEX (mm)	L (mm)	CODE
	1.25mm	7	HD1207
		10	HD1210
		15	HD1215
		25	HD1225

FOR SLIM IMPLANT

	HEX (mm)	L (mm)	CODE
	2.0mm	10	HD2010-SL
		15	HD2015-SL

FOR INTERNAL HEX IMPLANT

	HEX (mm)	L (mm)	CODE
	2.42mm	7	HD2407
		10	HD2410
		15	HD2415
		25	HD2425

FOR CONICAL NP IMPLANT



	HEX (mm)	L (mm)	CODE
	2.25mm	10	HD2210
		15	HD2215

FOR CONICAL RP IMPLANT


	HEX (mm)	L (mm)	CODE
	2.65mm	10	HD2610
		15	HD2615

FOR FLEXIBLE NECK ONE PIECE IMPLANT

	L (mm)	CODE
	18	RDOP



MOTOR MOUNTS



FOR SLIM EVA MOUNTLESS IMPLANT

	HEX (mm)	L (mm)	CODE
	2.0mm	15	HD2015-ML

FOR EVA MOUNTLESS IMPLANT

	HEX (mm)	L (mm)	CODE
	2.42	10	HD2410-ML
		15	HD2415-ML

FOR ONE PIECE IMPLANT

	L (mm)	CODE
	15	HDOP
	18	HDO
	50	HDOL

FOR MULTI UNIT 1.4

	L (mm)	CODE
	6	NM-HD06
	10	NM-HD10

HAND HEX DRIVER FOR MULTI UNIT 1.4

	L (mm)	CODE
	23.65	NM-HHD

FOR PROSTHETICS

	HEX (mm)	L (mm)	CODE
	1.25mm	20	MM1220
		28	MM1228

FOR SLIM IMPLANT

	HEX (mm)	L (mm)	CODE
	2.0mm	20	MM2020-SL
		28	MM2028-SL

FOR INTERNAL HEX IMPLANT

	HEX (mm)	L (mm)	CODE
	2.42mm	20	MM2420
		28	MM2428

FOR SLIM EVA MOUNTLESS IMPLANT



HEX (mm)	L (mm)	CODE
2.0	22	MMIP22-SL
	28	MMIP28-SL

FOR EVA MOUNTLESS IMPLANT



HEX (mm)	L (mm)	CODE
2.42	22	MM2422-ML
	28	MM2428-ML

FOR CONICAL NP IMPLANT



HEX (mm)	L (mm)	CODE
2.25mm	22	MM2220
	28	MM2228

FOR CONICAL RP IMPLANT



HEX (mm)	L (mm)	CODE
2.65mm	22	MM2620
	28	MM2628

FOR ANGULATED CAD/CAM



HEX (mm)	L (mm)	CODE
Ø1.7	24	MMAT

UNIVERSAL MOTOR MOUNT FOR IMPLANTS & PROSTHETICS


FOR INTERNAL HEX IMPLANT




HEX (mm)	L (mm)	CODE
2.42	25	MMIP20
	31	MMIP24

UNIVERSAL DRIVER KIT

FOR INTERNAL HEX IMPLANT



TYPE	L (mm)	CODE
2.42	31	UNIP24
Hex Driver Adapter	22	
Manual Driver Adapter	15.5	



LONG HANDLE DRIVERS

FOR PROSTHETICS



HEX (mm)	L (mm)	CODE
1.25mm	60	HHD1260

FOR INTERNAL HEX IMPLANT



HEX (mm)	L (mm)	CODE
2.42mm	143	HD

FOR ONE PIECE IMPLANTS



L (mm)	CODE
175	HDLOF

MANUAL WRENCH DRIVER ADAPTER



HEX (mm)	CODE
6.35	HRAL

MANUAL TORQUE DRIVER

For tightening of Cover Screws and Healing Caps.

MATERIAL: Stainless Steel. **TORQUE LIMITATION:** 15 Ncm.

LONG HANDLE



HEX (mm)	CODE
1.25	MTDH

SHORT HANDLE



HEX (mm)	CODE
1.25	MTD

GDT ELECTRIC ANGLE DRIVER

Combines a powerful motor with an optimized 90° angle head, ensuring the perfect torque for precise abutment fixation, even in tight spaces.

TORQUE OPTIONS: 20 Ncm, and 30 Ncm.

CHARGING TIME: 50 Minutes.

NET WEIGHT: 74g.

SCREWDRIVER WEIGHT: 56g.

PACKAGED WITH: Electric Driver, Short Bit - 1.25mm, Long Bit - 1.25mm, Universal Adapter, C-Type Charging Cable, Instructions.



NAME	CODE: EADP	
Short Bit		1.25mm
Long Bit		
Universal Adapter		

INSTRUMENTS

TOOLS

Wrenches

108

Adapters

108

Laboratory


108

Osteotomes

108

Measuring Tools

109



STAINLESS STEEL

ADAPTER

Ø6.35mm

WRENCHES

RATCHET WRENCH	HEX (mm)	CODE
	6.35	RA
TORQUE RATCHET WRENCH 10-45 Ncm	HEX (mm)	CODE
	6.35	TRA
TORQUE RATCHET WRENCH 25-100 Ncm	HEX (mm)	CODE
	6.35	STRA
HAND WRENCH	HEX (mm)	CODE
	6.35	WH

ADAPTERS

MOTOR MOUNT ADAPTER FOR RATCHET	HEX (mm)	CODE
	6.35	MMRB
HEX DRIVER ADAPTER FOR MOTOR MOUNT	HEX (mm)	CODE
	6.35	RDMM

LABORATORY

MANUAL ANALOG DRIVER	HEX (mm)	CODE
	2.42	AHD

OSTEOTOMES

PTERYGOID BONE OSTEOTOMES	D (mm)	CODE
	Ø2	M-X200
	Ø3	M-X300

MEASURING TOOLS

TYPE	PARALLEL PIN	L (mm)	CODE
Short		10	PP10
Long		16	PP16
Double Sided		6 / 16	PP30

IMPLANT DEPTH PROBE	CODE
 Osteotomy Depth (mm) Gingival Height (mm)	DEGG

IMPLANT DEPTH PROBE WITH CASE	CODE
	BTIP0001

ZYGOMATIC IMPLANT DEPTH PROBE	CODE
	DEZG

ANGULATION GUIDE GAUGE CLIPS			 2.42
ANGLE	17°	30°	45°
			
CODE	AGP17	AGP30	AGP45

ALL SOLUTION FOR A PERFECT **S**mile



www.gdentalt.com

Global Dental Transfers Ltd.

M Tower Building
Recha Freier 9
Beer Sheva
8489440
Israel

Tel: +1 929 442 0773
Email: Info@gdentalt.com